# EUROMEDEX

## Carbonate-bicarbonate Coating Buffer Tablets pH 9.6

Article no.	Product name	Pack size	Solution vol.
EU1-8932-50	Carbonate-bicarbonate buffer, pH9.6 with azide	50 tablets	100 ml/tablet
EU1-9832-8	Carbonate-bicarbonate buffer, pH 9.6 with azide	8 tablets*	100 ml/tablet
EU1-8922-100	Carbonate-bicarbonate buffer, pH 9.6	100 tablets	100 ml/tablet
EU1-8922-24	Carbonate-bicarbonate buffer, pH 9.6	24 tablets*	100 ml/tablet
Blister pack			

### Features

- Supplied with or without sodium azide as preservative
- Dissolve and use in minutes
- Exactly pre-weighed tablets
- Increased precision
- Pre-set pH

### **Product description**

Carbonate bi-carbonate buffer is used as coating buffer on plates or tubes. Carbonate bi-carbonate buffer is specifically developed for protein-coating procedures on microtiter plates and plastic tubes for RIA and EIA techniques. The buffer is supplied in two formulations, one of which contains 0.05 % sodium azide as preservative.

The buffer is supplied as exactly pre-weighed tablets giving 100 ml of 0.05 M sodium carbonate-bicarbonate buffer with pH 9.6 at 25 °C when dissolved in deionized water.

### Applications

• Protein coating on microtiter plates and plastic tubes.

### Directions for use

Deposit one tablet in a laboratory flask or beaker placed on a magnetic stirrer. Add 50 ml of deionized water and stir for a couple of minutes. Adjust the volume up to 100 ml, stir until full dissolution and the buffer is ready to use.

### Tips and hints

- If tablet is not properly dissolved, make sure:
- the water temperature is 22-25°C
- the buffer is properly stirred.

# Specifications Chemicals Analytical grade Format Exactly pre-weighed tablets Concentration 0.050 M sodium carbonate bicarbonate buffer, pH 9.6 at 25°C Volume 100 ml pH 9.6 ± 0.05 at 25°C Shelf life Three years after production date

### Shipping and storage

Carbonate-bicarbonate buffer is shipped at room temperature. Store tablets in a dry place at room temperature. Shelf life is three years.

### Certifications

Each stage of the manufacturing process is controlled and monitored by stringent quality control procedures to guarantee the highest possible quality and lot-to-lot reproducibility.



24, rue des Tuileries BP684 67460 SOUFFELWEYERSHEIM Tél. : 03 88 18 07 22 Fax.: 03 88 18 07 25 www.euromedex.com / research@euromedex.com